

GD/HS19 Diabetic Foot Nursing Simulator

GD/HS19



Features:


1. Designed for nursing of foot pathological changes caused by diabetes
2. No. 1,2,3 invasive trauma with minor foot infection
3. Display severe foot diseases, such as toe amputation, charcot foot and gangrene etc.
4. Soft and elastic model material, bendable toes

GD/L63 Prostate Inspection Model

This model provides prostates with tumors in 4 different development periods.

1. The 1st development :Beginning, a little big but the prostates still works well.
2. The 2nd development : Some hard tubers can be bserveed in the top right quadrant of prostate. It's the beginning of the tumors.
3. The 3rd development : The tumor is proliferating, and the previous tubers become bumps.
4. The 4th development : The prostates is all replaced by tumors, the surface become more rough and harder than before.

GD/L63



GD/L64 Peritoneal Dialysis Simulator

Function :

1. This model is a specially designed simulator for the medical staff to practise the skill of peritoneal dialysis.
2. Fistula orifice is in the realistic position, and the simulative liquid flows out through it.

Standard Equipment :

- a. Catheter
- b. Disposable Water-Proof Fabric
- c. Infusion Stand
- d. Liquid Bag
- e. Deluxe Portable Box

GD/L64



GD/LV1 Advanced Surgical Suture Arm

This model of Advanced Surgical Suturing Arm allows the students to practise cut, suturing, stitches release, enswathing, etc. The elasticity and flexibility of skin can make it for practise hundreds of times. There are some initial wounds on the arm if all the initial wounds have been sutured, people can cut the arm in other place for new training of suturing.

GD/LV1

