

GD/A11301/1 Full Size Human Body Showing Muscles & Organs (86 Cm) :

This full size high quality muscular figure represents a wide variety of human anatomical structures with accurate details. Hence its suitable for various medical colleges, universities, nursing training institutes, laboratories & hospitals.


This model consists of 27 parts, such as muscles of chest wall and abdomen, muscles of upper and lower limbs, skull, brain and viscus. It shows the structure of head, neck, torso, upper and lower limbs, muscles, muscular tendon, ligaments, viscus, blood vessels, brain, etc.

Material : Advanced PVC and painted with imported paint.




MK-050 Human Torso

- ◆ Hand painted in vivid and natural colours
- ◆ 11 dissectible parts
- ◆ ABS Base
- ◆ Fully illustrated colour
- ◆ Instructional manual



MK-064 56cm (22") visible Expectant Mother Anatomy Kit :

- ◆ Complete Skin, Skeleton and Vital Organs
- ◆ Learn the different part of the Human Body
- ◆ Easy Assembly
- ◆ Illustrated Instruction Manual



XC-204 Unisex Torso 23 Parts (85 cm)

This outstanding torso features an exposed spine with removable vertebrae and spinal cord segments, a female breast plate and interchangeable male and female genitalia. The female organs include a fetus in the womb. Dissectible into 23 parts : torso, female breast plate, head, eyeball, brain, vertebra, spinal nerves, lung (2 parts), heart (2 parts), liver, kidney, stomach (2 parts), intestines (4 parts), male genitalia (2 parts), female genitalia with fetus (3 parts).

Made of PVC plastic. Mounted on a plastic base.

